

## Recommend Spare Parts for NMR-271

	<b>Description</b>	<b>Part No.</b>	<b>Field stock</b>	<b>Service Co.</b>
1	Boot, SW1	K24-17239-00	2	5
2	Boot, Circuit Breaker	K24-18164-00	2	5
3	Circuit Breaker, 25A (CB1)	K24-22363-25	1	2
4	Circuit Breaker, 8A (CB2)	K24-20565-00	1	2
5	Circuit Breaker, 4 Amp (CB3)	K24-22330-00	1	2
6	Switch, Manual Defrost (SW2)	K24-21355-00	1	2
7	Switch, On-Off (SW1)	K24-22144-00	1	2
8	Contactora Compressor	k24 22403 00	1	2
9	Contactora Reversing	k24 22169 00	0	2
10	Contactora Condensor-Evaporator	k24 22155 00	1	2
11	Transformer 460 – 12/24/110V	k25 26404 00	1	2
12	Coil, Solenoid Valve	K25-26103-01	1	2
13	Thermostat	K31-00887-10	1	3
14	Phase Sensor	K31-00809-00	1	2
15	Touch Pad	K31-00811-00	1	2
16	Switch, High Pressure	360-13230-01	1	2
17	Ribbon Cable 18" long	360-12540-01	1	2
18	Battery Charger	360 12473 11	1	2
19	Probe, Kit return & Defrost NMR	360 13807 01	1	2
20	Switch, Discharge Temperature	360 13234 01	1	2
21	Evaporator Fan Motor	K24-22162-00	0	2
22	Condenser Fan Motor	K24-22316-00	1	2
23	Compressor	K15-00043-00	0	2